# VERI-Q MULTI Monitoring System

**Sample Volume:** 15µL for Lipids 7µL for hemoglobin 0.5 µL for glucose

### **Test Speed:**

3 minutes for Lipids 5 seconds for glucose and hemoglobin

#### Advantages:

- The Veri-Q needs only 15 µL of whole blood for the enzymatic reaction to occur, providing a more precise and accurate result in comparison with other POC devices, and automatic or semi-automatic analyzers
- Rapid and accurate results with low test volume
- Stores up to 1,000 results
- Reading method: electro-chemical for glucose and photometric for lipids and hemoglobin

## Features

Extremely accurate and simple whole blood cholesterol analyzer, not bigger than a mobile phone. Test your full lipid panel in less than 3 minutes; a super tool for monitoring HDL and LDL. Quantitative measurement of six parameters (TC/HDL/TRIG/LDL) + (TC & HDL Ratio) + (LDL & HDL Ratio)



Lipid Panel: Simultaneous measurement of six parameters (TC/HDL/TRIG/LDL) + (TC & HDL Ratio) + (LDL & HDL Ratio) Includes: Memo Chip and Micropipettes for blood collection. The lipid tests are also supplied separately for specific tests.

Glucose test strips for the Veri-Q Multi Meter Analyzer. **Pack Size**:50 Strips





600 W 21 Street, Hialeah, Florida 33010 USA Tel. (305) 826-9886 Fax (305)826-9887 WWW.IDGone.com